- de** Gebrauchsanweisung
- en** User's Manual
- fr** Mode d'emploi
- it** Istruzioni per l'uso
- es** Instrucciones de empleo
- pt** Manual de operação
- ru** Руководство по применению

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Dermatoskope

Kontaktscheibe mit Skalierung für Dermatoskope Ø 25 mm
 Contact disk with scale setting for dermatoscopes Ø 25 mm
 Rondelle de contact avec cadran pour dermatoscopes Ø 25 mm
 Vetrino di contatto con scala graduata per dermatoscopi Ø 25 mm
 Arandela de contacto con escala para dermatoscopios Ø 25 mm
 Disco de contacto com escala para dermatoscópios Ø 25 mm
 Контактное стекло со шкалой Ø 25 мм

Fokuseinstellung
 Focus adjustment
 Plage de focalisation
 Messa a fuoco
 Margen de la focalización
 Ajuste do foco
 Установка фокуса

Ladebatterie (Akku)
 Rechargeable battery
 Pile rechargeable (accu)
 Batteria ricaricabile
 Pila recargable (acumulador)
 Pilha recarregável (acumulador)
 Аккумулятор

Batterien
 Batteries
 Piles
 Batterie
 Pilas
 Pilhas
 Батарейки

de	Dermatoskop	4
en	Dermatoscope	9
fr	Dermatoscope	13
it	Dermatoscopio	18
es	Dermatoscopio	23
pt	Dermatoscópio	28
ru	Дерматоскопы	33

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Dear buyer! Thank you for choosing a KaWe product. Our products are known for high quality and durability. This KaWe device complies with the provisions of Regulation EU 2017/745 (European Regulation on Medical Devices), whereby it is classified as a medical product of class 1.

Please carefully consider this instruction and read all the provisions prior to use. Follow the maintenance recommendations.

Prior to use, please carefully consider the maintenance procedure.

Application: Only qualified medical professionals can use the dermatoscope.

Group of patients with indications for use:

People of any age, male and female, requiring an examination of skin surface structures.

Purpose: KaWe dermatoscope is designed for visual examination of body surface. It zooms skin surface structures almost ten-fold, and with a halogen lamp lights the area in question. The mechanism is analogous to microscopic examination in reflected light (epiluminescence microscopy). Dermatoscopy is aimed at precise examination of skin pigment changes in everyday practice, specifically, within early detection of malignant melanoma. Skin reflexes are suppressed by a special dermatoscopic moisturizing gel (KaWe Derma-Gel), skin color changes are localized and detected more efficiently. Contact lens scale allows establishing the sizes of changes.

The device is designed only for skin surface examination.

Dermatoskope:

- **KaWe EUROLIGHT** D30 2.5 V dermatoscope
- **KaWe PICCOLIGHT** D 2.5 V dermatoscope
- **KaWe EUROLIGHT** D30 LED 3.5 V dermatoscope

Maintenance of KaWe EUROLIGHT D30 2.5 V and 3.5 V dermatoscopes:

Proper functionality is guaranteed when using KaWe standard handles with KaWe rechargeable batteries 2.5 V (for handles with 2.5 V bulbs) or 3.5 V (for handles with 3.5 V bulbs). For 2.5 V handles, it is also possible to use 2 standard batteries (type C / LR14) instead of rechargeable batteries. Put them on and fix the device's head. To switch the device on, press and turn to the left a white plastic button on the rotary switch. Variable resistor allows quick brightness regulation (0=off, 1=max. brightness). To switch the device off, turn the button the other way back to reference position (switching off movement). To replace the battery or rechargeable battery, open the lower cover. To charge handles with rechargeable batteries, take off the plastic cover and connect with KaWe Med-Charge 4000 charging device.

To charge KaWe EUROLIGHT D30 LED 3.5 V, unscrew two connecting parts of the handle and put the lower part of the handle with the mains plug into the socket. (ATTENTION: 3.5 V version is compatible with KaWe 3.5 V rechargeable batteries only).

Remove the investigation head by lowering the throw-on insert.

ATTENTION: Follow additional maintenance recommendations to make the service life of the rechargeable batteries longer.

Maintenance of KaWe PICCOLIGHT: Proper functioning is provided by KaWe PICCOLIGHT handles with 2 batteries (Mignon type AA / LR6). To switch the device on/off, use a sliding switch. To switch the device on using a sliding switch, move it towards the dermatoscope's head. To switch the device off, move the sliding switch in reverse direction, to reference position. To replace the battery and lamp, unscrew the head of KaWe PICCOLIGHT.

Time of application: The devices are designed for periodic use, maximum operation time totals 1 minute. There shall be 10-minute pauses between the applications.

Lamp replacement of KaWe PICCOLIGHT: For proper light functioning and quality, use only original KaWe lamps. Unscrew the head of the dermatoscope from the handle and remove the lamp. Clean the glass bulb of the new lamp with alcohol-containing cloth. The glass bulb shall be clean and free from fingerprints (grease spots). Put the lamp all the way into its groove in the side pin. Then screw the investigation head back to the handle.

Lamp replacement of KaWe EUROLIGHT: For proper light functioning and quality, use only original KaWe lamps. Unlatch the head of the dermatoscope from the handle and remove the lamp. Now repeat the steps above in the section on replacing KaWe PICCOLIGHT lamp. The glass bulb shall be clean and free from fingerprints (grease spots). Put the lamp all the way in.

Handling: Put special moisturizing KaWe Derma-Gel dermatoscope gel on the skin area of interest. Turn the light on and set the contact lens on skin without excessive pressing. Even lighting will cover the skin area with large examination diameter (25 mm). Place the examined area in the center of the contact lens, if possible. Bring the eyepiece as close to the eye as possible. With a free hand, turn the focus ring until the image of the area of interest becomes sharp. Try looking through the eyepiece without pressure on eye muscles (not accommodating them for sharper view). Scales contact lens is designed to control the pigmented skin changes. To examine raised or palpable skin areas (damages), dermatoscope can be used without the contact lens. For a sharper view you can raise the dermatoscope over the area of interest but not further than 5 mm.

Technical characteristics:

- Magnifier glass with 10x magnification
- Contact lens diameter: 25 mm
- Focusing / sharpness adjustment: -6...+3.5 diopter
- Halogen lamp aperture: 2.5 V
(KaWe EUROLIGHT D30)
- Vacuum lamp aperture: 2.5 V
(KaWe PICCOLIGHT D)
- LED lamp aperture: 3.5 V
(KaWe EUROLIGHT D30 LED)

Additional directions, maintenance and storage: The device will serve for years with proper handling and storage.

Warranty: When used properly and with attention to our user's manual as well as when used with original parts, we warranty this product for two years from the date of purchase. If you have any questions, or your device requires repair, please contact your supplier.

Warnings: An error may cause extensive heat posing a threat of burn injury. If the device heats excessively, immediately turn the handle off.

Warnings: Do not put the dermatoscope in fluid and keep fluid off the handle! The dermatoscope is not designed for surgery, so it is sufficient to clean or disinfect it with surface alcohol-containing disinfectants. Do not use abrasive cleaning agents! The dermatoscope shall not be sterilized.

Limitations of reprocessing: Reprocessing has a minor effect on these devices. The service life of the product is usually defined based on the degree of wear and use-related damage.

Directives

Working area: Clean the surface with a disposable/paper cloth.

Environment conditions:

	TEMPERATURE	RELATIVE AIR HUMIDITY
OPERATION	+10°C to +35°C	max. 80 %
STORAGE	-25°C to +70°C	max. 95 %
TRANSPORT	-25°C to +70°C	max. 95 %

Transportation and storage: Store and transport dermatoscopes protecting them from dust, moisture and dirt.

Preparation for cleaning: No special requirements.

Manual cleaning: Clean the dermatoscope external surface with a damp, soft and lint-free cloth. Do not put the dermatoscope in fluid and keep body off the handle! The dermatoscope is not designed for surgery, so it is sufficient to clean it with surface alcohol-containing detergent. Do not use abrasive cleaning agents! The device shall not be sterilized. Clean the magnifying glass with a glass detergent only.

Automatic cleaning: Automatic cleaning is impossible/unnecessary.

Disinfection: Clean the dermatoscope with surface alcohol-containing disinfectants. Use a slightly damp cloth.

Sterilization: Contact disk shall be sterilized in autoclave at 134 °C for at least 3 minutes.

Additional information: When sterilizing several devices in autoclave, do not exceed the maximum load indicated by the sterilizer manufacturer.

Maintenance: Remove or replace the damaged parts. If the device does not work, do the following: **1.** Replace the batteries or rechargeable batteries and check whether they are serviceable; **2.** Replace the lamp and check whether it is serviceable; **3.** Send the device for repair.

Control and functional check: Check the smooth travel of moveable parts, such as focus adjustment frame and plugs. Switch the device on. All the instruments: Visually check for damages and wear. Keep all the electric contacts clean.

Package: Separate zipped bag / fabric bag.

Disposal: Electric/electronic devices shall not be disposed together with household waste! Dispose the device in accordance with the local regulations at a special collection point.

General accessories: To find out more about these products, visit our homepage: www.kawemed.com.

Manufacturer: **KaWe**

Manufacturer's contact details: Dealer's address or phone number: +49-7141-68188-0.

Basic UDI-DI

4030155KaWe0103YV

Legend:

	Manufacturer & date of manufacture
	Heed user's manual
	Safety instructions or warning
	Separate collection of electrical / electronic equipment and rechargeable batteries
	Item No.
	CE conformity label
	Applied part Type BF
	Store in dry conditions
	Temperature range limits
	Medical product
	UDI Data Carrier
	Use by date

Instructions for service personnel and patients

In case of problems during the device use, it is necessary to immediately contact the manufacturer and the competent authorities of the relevant member state of the Union, in which the service personnel or the patient are located.

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KIRCHNER & WILHELM GmbH + Co. KG