
de	Gebrauchsanweisung	MASTERLIGHT (Untersuchungsleuchte)
en	User's Manual	MASTERLIGHT (examination light)
fr	Mode d'emploi	MASTERLIGHT (lampe d'examen)
it	Istruzioni per l'uso	MASTERLIGHT (lampada da osservazione)
es	Instrucciones de empleo	MASTERLIGHT (lámpara de exploración)
pt	Manual de operação	MASTERLIGHT (lâmpada de observação)
ru	Руководство по применению	MASTERLIGHT (смотровой светильник)

de

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en

Please read this User's Manual thoroughly and carefully before attempting to use this product and heed the given care instructions.

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de Montage MASTERLIGHT LED 1000

en Installation of MASTERLIGHT LED 1000


fr Installation de MASTERLIGHT LED 1000

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





















5. Insert the lamp shade. Now tighten the screw (9) to prevent accidental removal of the lamp shade.

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
6. Height adjustment: using your thumb, press the black locking ring (10) upwards and set the desired height.

Attention: Lowering operation stops immediately after releasing the ring!

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Table of Contents

1. Symbols used in the User's Manual and on packaging	29
2. Purpose and intended use	29
3.  Safety instructions / contraindications	29
4. Installation and maintenance	30
MASTERLIGHT LED 1000	30
A) Scope of delivery	30
B) Installation	30
C) Maintenance: ON/OFF switch	30
D) Maintenance: Positioning	30
MASTERLIGHT LED 2000	30
A) Scope of delivery	30
B) Installation	30
C) Maintenance: ON/OFF switch	31
D) Maintenance: Brightness control	31
E) Maintenance: Positioning	31
MASTERLIGHT LED 3000 / 3000F	31
A) Scope of delivery	31
B) Installation	31
C) Maintenance: ON/OFF switch	31
D) Maintenance: Brightness control	32
E) Maintenance: Focusing (3000F only)	32
F) Maintenance: Positioning	32
MASTERLIGHT LED 4000 / 4000F	32
A) Scope of delivery	32
B) Installation	32
C) Maintenance: ON/OFF switch	32
D) Maintenance: Brightness control	32
E) Maintenance: Focusing (4000F only)	32
F) Maintenance: Positioning	32
5. Preparatory measures before use	32
6. Hygienic preparation	32
Manual cleaning/disinfection	32
A) Stand	33
B) Lamp shader	33
C) Protective glass	33

D) Handle (except for MASTERLIGHT LED 1000)	33
Automatic cleaning.....	33
Sterilization	33

7. Storage and transportation	33
Ambient temperature.....	33
8. Maintenance and service	33
Maintenance.....	33
Clamping lever adjustment	34
9. Technical data	34
Electrical installation information.....	34
Other important information	35
10. Accessories and spare parts	35
11. Disposal	35
12. Appendix	35
Warranty	35
Contact information.....	35
Guidelines and tables on electromagnetic compatibility	35

Dear Customer, thank you for choosing a KaWe product. Our examination lights feature powerful LEDs. They comply with the provisions of Regulation (EU) 2017/745 (European Medical Device Regulation) and, according to this Regulation, are classified as Class I medical products.

 Please read and follow this User's Manual, and keep it in case you need to read it again. Before use, thoroughly read the maintenance procedure and follow the care instructions.

In case of uncertainty about the connection or maintenance of the device, please contact the customer service (see the manufacturer's address on the back).

KaWe offers the following examination light models:

- MASTERLIGHT LED 1000
 - MASTERLIGHT LED 2000
 - MASTERLIGHT LED 3000
 - MASTERLIGHT LED 3000F (focusable)
 - MASTERLIGHT LED 4000
 - MASTERLIGHT LED 4000F (focusable)
- (See Technical Data in Section 9)

The advantages of LED technology:

- Long service life of min. 60.000 hours
- Low power consumption
- Heat radiation reduced to a minimum
- Natural and contrast color rendering
- Accurate highlighting of the damaged area

1. Symbols used in the User's Manual and on packaging

	Caution
	Please follow the User's Manual (color: blue)
	Product reference number
	Series number
	Manufacturer & if applicable date of manufacture
	Date of manufacture
	Protection class I
	Protection class II

	Keep dry
	Fragile, requires careful handling
	Temperature range limits
	Humidity, limit
	This side up
	Separate collection of electrical and electronic equipment
	Non-ionizing radiation
	CE mark of conformity
	Medical device
	UDI Data Carrier

2. Purpose and intended use

The MASTERLIGHT examination lights are designed to aid treatment and diagnosis in a hospital or medical office. They illuminate the work area and the patient's body with a shadow-free, "cold" high power light.

- Patient group for which use is indicated: all patient groups
- This product should only be used by authorized trained personnel!

Compliance with all instructions in this User's Manual as well as important regulations and standards also applies to the intended use. Any other use is not intended; therefore, the manufacturer is not liable for any damage resulting from such use!

3. Safety instructions / contraindications

- **Be careful when taking the stand out of the box!** The stand is live so there is risk of electric shock. The inner tube is extendable.
- This device is not designed for use in explosion hazardous areas!
- These devices are examination lights in accordance with EN 60601-2-41, which are not fail-safe free-standing lighting devices.
- Protect the examination light from moisture!
- The light must only be mounted on a stand!
- Before installation, leave the device in the packaging for at least 24 hours in a suitable room to compensate for temperature changes.

- The socket used must be installed in accordance with the requirements of IEC or VDE 0100-710.
- Examination light repairs, and especially installation work, must only be carried out by KaWe or by an organization specifically authorized by KaWe.
- KaWe is only responsible for the safety of the stand and examination light when repairs and modifications are made by KaWe itself or by an organization ensuring that the safety regulations are observed.

Modifications to examination light are **prohibited!**

- Use only original KaWe accessories and spare parts.
- The manufacturer will not be held liable for damage to health or property if the stand or examination light is not maintained properly, not maintained in accordance with their purpose, or misused.
- The examination light must be maintained every two years.
- The connection to the mains is made via a plug cable.

4. Installation and maintenance

MASTERLIGHT LED 1000

A) Scope of delivery

- 1x stand base (five-footed base with foot juncture piece, foot pieces including rollers and hardware)
- 1x stand post
- 1x lamp shade incl. arm
- 1x User's Manual


B) Installation

See "Installation of MASTERLIGHT LED 1000" on page 4.

During the installation of the examination light, the entire system must be disconnected from the mains!

C) Maintenance: ON/OFF switch


To turn the examination light on and off, press the (A) switch on the lamp shade.

D) Maintenance: Positioning


Press the clamping lever (K) on the arm.

The lamp shade can be moved in the indicated directions. The flexible swivel (G) holds the lamp shade in position.

The examination light must not be installed at an angle $< 90^\circ$.

Risk of rupture of flexible hose connection!

MASTERLIGHT LED 2000

A) Scope of delivery

- 1x stand base (five-footed base with rollers, protective cap and hardware)
- 1x stand post
- 1x lamp shade incl. arm
- 1x User's Manual


B) Installation

See "Installation of MASTERLIGHT LED 2000" on page 7.

During the installation of the examination light, the entire system must be disconnected from the mains!

C) Maintenance: ON/OFF switch

To turn the examination light on and off, press the (A) switch on the power supply housing.

D) Maintenance: Brightness control

Brightness is controlled on the control panel (B) in the range from 50% to 100%:

- ⇒ Press the „Intensity +„ button: Brightness will increase
- ⇒ Press the „Intensity -„ button: Brightness will decrease

E) Maintenance: Positioning

Press the clamping lever (K) on the transformer housing. Use the handle (H) to adjust the lamp shade position. The flexible swivel (G) holds the lamp shade in position.

The examination light must not be installed at an angle $< 90^\circ$.

Risk of rupture of flexible hose connection!

MASTERLIGHT LED 3000 / 3000F

A) Scope of delivery

- 1x stand base (five-footed base with foot, rollers and hardware)
- 1x stand post
- 1x lamp shade incl. short arm
- 1x User's Manual


B) Installation

See Installation of MASTERLIGHT LED 3000/ 3000F & MASTERLIGHT LED 4000/ 4000F on page 10.

During the installation of the examination light, the entire system must be disconnected from the mains!

C) Maintenance: ON/OFF switch


To turn the examination light on and off, press the (A) button on the control panel.

To switch off, press and hold button (A) for at least 1 sec. If the light is in standby mode, the top LED on the control panel flashes.

D) Maintenance: Brightness control

Brightness is controlled on the control panel in the range from 50% to 100%:

- ⇒ Press the „+“ button (D): Brightness will increase
- ⇒ Press the „-“ button (B): Brightness will decrease


The set intensity is displayed on the (C) indicator.

E) Maintenance: Focusing (3000F only)

To focus (change the light spot size) turn the handle (H) on the lamp shade (see figure).

F) Maintenance: Positioning

The lamp shade position can be adjusted with the handle (H) or the two handle sleeves (G).

Before each use make sure the examination light is in good working order!

MASTERLIGHT LED 4000 /4000F

A) Scope of delivery

- 1x stand base (five-footed base with foot, rollers and hardware)
- 1x stand post
- 1x lamp shade with short arm
- 1x User's Manual

B) Installation

See Installation of MASTERLIGHT LED 3000/ 3000F & MASTERLIGHT LED 4000/ 4000F on page 10.

During the installation of the examination light, the entire system must be disconnected from the mains!

C) Maintenance: ON/OFF switch

See MASTERLIGHT LED 3000/3000F on page 31.

D) Maintenance: Brightness control

See MASTERLIGHT LED 3000/3000F on page 32.

E) Maintenance: Focusing (4000F only)

See MASTERLIGHT LED 3000/3000F on page 32.

F) Maintenance: Positioning

See MASTERLIGHT LED 3000/3000F on page 32.

5. Preparatory measures before use

During initial commissioning, the device must be checked in accordance with EN 62353!

In particular, check the examination light for the following points:

- Paint defects
- Cracks in plastic parts
- Loose parts
- Check the connection between the examination light and support system
- Fail-safe operation (e.g., brake rollers)
- Electrical safety

Defective parts must be replaced before continuing to use. Tighten the screws if necessary.

6. Hygienic preparation

Clean and disinfect immediately after use. Remove dirt from the rollers to prevent static charge build-up.

Manual cleaning / disinfection

For disinfection, it is recommended to use products dissolved in water or in max. 20% alcohol.

Do not use acid-containing products!

**Protect the device from moisture penetration!**

Do not allow liquid to get inside the housing.

A) Stand

The surface of the stands is easily maintained by damp wiping. Common disinfectants can be used for cleaning.

B) Lamp shade

The examination light has a high quality surface. The cleanliness of the light surface is easily maintained by damp wiping with a common disinfectant.

C) Protective glass

The protective glass is made of high quality polymer material. It is possible to use common disinfectants to clean it. Observe the following instructions for cleaning:

- ⇒ Always wipe the protective glass (S) with a damp material (never with dry)!
- ⇒ After cleaning, wipe the protective glass (S) with an antistatic agent. Use a lint-free cloth for this.

**D) Handle (except for MASTERLIGHT LED 1000)**

Common disinfectants can be used to clean the handle.

We recommend replacing the handle with a new one at least every two years. This ensures sustainable, stable and safe operation with the handle as well as cleaning/disinfection efficiency.

Automatic cleaning

Not applicable.

Sterilization

Not applicable.

7. Storage and transportation

During storage and transportation, examination lights must be protected from dust, moisture and bacterial contamination.

Environmental conditions

	Temperature	Relative humidity	Air pressure
Operation	from +10°C to +30°C	from 30% to 75%	from 700 hPa to 1060 hPa
Storage	from -10°C to +50°C	from 20% to 90%	
Transportation			

**8. Maintenance and service**

The end of product service life is usually determined by the degree of wear and tear due to use. Damaged parts must be rejected and replaced!

For repairs, please contact your dealer.

Maintenance

Maintenance and inspection of the examination light must be carried out **at least every 2 years**. This includes technical and mechanical inspection.

The inspection of the light and the supporting structure is carried out on the basis of DGUV V3 (formerly BGV-A3) in connection with EN 62353.

Before any maintenance or inspection operations, switch off the examination light and unplug it! Prevent the light re-switching.


Clamping lever adjustment

MASTERLIGHT LED 1000 & MASTERLIGHT LED 2000

If the lamp shade moves too tightly or does not maintain the set position, it is necessary to adjust the braking effect of the clamping arm between the light arm and the stand post.

MASTERLIGHT LED 1000

MASTERLIGHT LED 2000

9. Technical data

Note: The technical data are subject to certain deviations. For operational reasons, the actual values may differ slightly from the values specified below. The values for Ra can deviate by approx. $\pm 5\%$. The color temperature values can deviate by APPROX. ± 200 K.

	MASTERLIGHT model					
	LED 1000	LED 2000	LED 3000	LED 3000F	LED 4000	LED 4000F
Photometric data						
Brightness / 0.5 meter [lux]	30 000	60 000	-	-	-	-
Brightness / 1 meter [lux]	7 800	14 000	40 000	50 000	60 000	70 000
Color rendering index Ra	97	97	96	96	96	96
Color rendering index R9	96	95	94	94	94	94
Color temperature [K]	4000	4000	4500	4500	4500	4500
Heat build-up (head) [°C]	0.5	0.5	0.5	0.5	0.5	0.5
Light spot size Ø [cm]	10	11	15	13-18	14	13-19
Dimming range [%]		50-100	50-100	50-100	50-100	50-100
Number of LEDs	3	7	12	12	19	19
Working lifespan of the LEDs [h]	60 000	60 000	60 000	60 000	60 000	60 000
Work area [cm]	35-80	40-100	70-140	70-140	70-140	70-140
Lamp shade Ø [cm]	12	22	29	29	33	33
Height adjustment [cm]	60	60	121	121	123	123
Electrical data						
Power input [W]	7	10	9	9	16	16
Operating voltage [V/DC]	12	24	24	24	24	24
Input voltage [V/AC]	100-240	90-264	100-240	100-240	100-240	100-240
Input frequency [Hz]	50-60	50-60	50-60	50-60	50-60	50-60
Amperage [A]	0.6	0.42	0.4	0.4	0.7	0.7
Protection class	II	I	I	I	I	I
Lamp shade protection class (acc. to IEC60529)	IP 42	IP 42	IP 42	IP 42	IP 42	IP 42
Weight						
Weight of lamp shade + stand [kg]	5	6	12	12	12.5	12.5

The rest of the technical data can be requested from the manufacturer (see the address on the back cover page).

Electrical installation information

When switched on, KaWe examination lights are exposed to a peak current. Their standard scope of delivery includes a power adapter with a plug.

The **MASTERLIGHT LED 2000/ 3000/ 3000F/ 4000/ 4000F** examination lights are devices with Protection Class I. To avoid the risk of electric shock, these devices may only be connected to a power adapter with a grounding wire!

Other important information

When several KaWe examination lights are used together, the total irradiation intensity due to the intersection of the light spots of these examination lights can exceed 1000 W/m². This leads to the risk of excessive heat build-up in the light spot. The intersection of the light spots of several KaWe examination lights may result in the UV limit (< 400 nm) exceeding 10 W/m.

Note: If you wish, it is possible to request a factory electrical safety test report. For this, we only need the series number of the KaWe examination light, which requires a report. If additional devices **or KaWe examination lights** are connected during installation, Section 16 of **EN 60601-1:2013** must be applied; may also need a compliance check!

10. Accessories and spare parts

Accessories and spare parts can be found in our general catalogue at <https://www.kawemed.com/download/broschueren/> or from your dealer.

Sterilizable handle models are available on request only.

11. Disposal

The examination light does not contain harmful substances. At the end of the product's service life, the examination light components must be properly disposed of. Observe carefully the requirements for careful separation of materials. The electrical printed circuit boards must be brought to an appropriate recycling organization. The light housing and other components of the examination light are disposed of according to specific material requirements.

12. Appendix

Warranty

With proper use and adherence to our User's Manual, as well as the use of original KaWe spare parts, the statutory warranty is valid for a period of two years from the date of sale (excluding light sources).

If you have any questions or need repairs, please contact your dealer.

Contact information

Contact your dealer or manufacturer (see back cover for address).

Basic UDI-DI

4030155KaWe10SK

Instructions for user and patient

All serious incidents involving the product must be immediately reported to the manufacturer and the competent authorities of the relevant Member State in which the user and/or patient is located.

Guidelines and tables on electromagnetic compatibility.

KaWe examination lights require special precautions related to electromagnetic compatibility (EMC) and must be installed in accordance with the EMC guidelines included in the accompanying documents. The performance of the examination lights can be affected by portable and mobile HF telecommunication devices.

The use of non-KaWe accessories or spare parts will result in increased emissions or reduced immunity of the device.

For proper operation of the examination light, it is necessary to prevent its placement in the immediate vicinity of other devices or on top of each other, and if it is required to operate it in the immediate vicinity of other devices or on top of each other, it is necessary to monitor the examination light.

Table 1: Environmental emission limits (according to IEC60601-1-2).

Phenomenon	Professional health care institutions	Home health environment
Conducted and radiated emissions	Class B, Group 1 (according to CISPR 11)	CISPR 11 ^{d,d)}
Distortion due to high-order harmonics	See IEC 61000-3-2 ^{b)}	See IEC 61000-3-2
Voltage fluctuations and flicker	See IEC 61000-3-3 ^{b)}	See IEC 61000-3-3
a) For information on categories of environments for intended use, see p. 8.9 b) These tests do not apply to this environment, unless the ME-DEVICES (medical electrical devices) and ME-SYSTEMS (medical electrical systems) used in it are connected to a public mains supply and their power consumption fluctuates in a range that is within the application range of the cited EMC basic standard. c) ME-DEVICES AND ME-SYSTEMS designed for use in aircraft must comply with the emission requirements according to ISO 7137. Measurement of conducted emissions is required only when the corresponding ME-DEVICES AND ME-SYSTEMS are designed to be connected to the aircraft on-board network. The ISO 7137 standard is identical to the RTCA DO-130C:1989 and EU-ROCAE ED-14C:1989 specifications. The latest editions of these specifications are RTCA DO-160G:2010 and EUROCAE ED-14G:2011. For this reason, Section 21 (and Category M) of a later edition, e.g., [39] or [40], should be considered. d) For other relevant operating conditions or other electromagnetic environments during transportation, the appropriate standards shall apply if they provide for the proper use of ME-DEVICES and ME-SYSTEMS. Examples of possible suitable standards include CISPR 25 and ISO 7637-2.		

Table 2: Shell (according to IEC60601-1-2).

Phenomenon	EMC basic standard or test method	Immunity test level	
		Professional health care institutions	Home health care environment
Static electricity discharge	IEC 61000-4-2	± 8 kV contact discharge ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air discharge	
High-frequency electromagnetic fields ^{a)}	IEC 61000-4-3	3 V/m ^{b)} 80 MHz to 2.7 GHz ^{b)} 80% AM at 1 kHz ^{c)}	10 V/m ^{b)} 80 MHz to 2.7 GHz ^{b)} 80% AM at 1 kHz ^{c)}
High-frequency electromagnetic fields in the vicinity of wireless communication devices	IEC 61000-4-3	See Table 9	
Magnetic fields with rated energy frequencies ^{d), e)}	IEC 61000-4-3	30 A/m ^{g)} 50 Hz or 60 Hz	

a) The interface between the patient physiological signal simulator, if used, and ME-DEVICES or ME-SYSTEMS should be located within 0.1 m from the vertical plane of the uniform field. The simulator must be directed in the same direction as the ME-DEVICE or ME-SYSTEM.

b) ME-DEVICES and ME-SYSTEMS, which for the purpose of their operation properly receive high-frequency electromagnetic energy, should also be tested at the receiving frequency. The test may also be conducted at other baseband frequencies as determined by the risk management process. This test evaluates the basic safety and essential performance of such a receiving ME-DEVICE in the presence of ambient signals at its receiving frequencies. It should be assumed that during this test, normal reception conditions for the receiver may not be achieved.

c) The test may also be conducted at other baseband frequencies as determined by the risk management process.

d) Applies only to ME-DEVICES with magnetically sensitive elements or connections.

e) When testing the ME-DEVICE, one of the provided rated input voltages can be applied, however, the corresponding mains frequency must correspond to the test signal frequency (see Table 1).

f) The level before modulation is applied.

g) This test level assumes the distance between the ME-DEVICE or ME-SYSTEM and the power frequency magnetic field sources of minimum 15 cm. If the risk analysis shows that ME-DEVICE or ME-SYSTEM, when used, turns out to be closer than 15 cm to the power frequency magnetic field sources, the test level of immunity should be adjusted in accordance with the expected minimum distance.

Table 3: AC 2-pole for power input (according to IEC 60601-1-2).

Phenomenon	EMC basic standard	Immunity test level	
		Professional health care institutions	Home health care environment
Electrical fast transient/ impulse withstand voltage a) l) o)	IEC 61000-4-4	± 2 kV 100 kHz repetition rate	
impulse voltages ^{a) b) j) k) o)} line-to-line	IEC 61000-4-5	± 0.5 kV, ± 1 kV	
impulse voltages ^{a) b) j) k) o)} line-to-earth	IEC 61000-4-5	± 0.5 kV, ± 1 kV, ± 2 kV	
Conducted interference caused by high frequency fields ^{c) d) e)}	IEC61000-4-6	3 V ^{m)} 0.15 MHz - 80 MHz 6 V ^{m)} in ISM bands from 0.15 MHz to 80 MHz ⁿ⁾ 80% AM at 1 kHz ^{e)}	3 V ^{m)} 0.15 MHz - 80 MHz 6 V ^{m)} in ISM and amateur radio bands from 0.15 MHz to 80 MHz ⁿ⁾ 80% AM at 1 kHz ^{e)}
Voltage dips ^{i) l)}	IEC 61000-4-11	0% U _T ; for 0.5 cycle ^{h)} at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° ^{o)} 0% U _T ; for 1 cycle and 70% U _T ; for 25/30 cycle ^{h)} Single phase: at 0 degrees	
Voltage interruption ^{f) i) o)} r)	IEC 61000-4-11	0% U _T ; for 25/30 cycles ^{h)}	
a) The test can be carried out at any mains input voltage in the range of RATED VALUES for ME-DEVICE or ME-SYSTEM mains voltage. If ME-DEVICE or ME-SYSTEM are tested at the mains input voltage, additional tests with other voltages are not required.			
b) During the testing, all the wires of ME-DEVICES and ME-SYSTEMS must be mounted.			
c) The calibration of the current supply clamps must be done on the 150.Ω system.			
d) If the frequency steps are outside the ISM or amateur radio band, an additional test frequency must be used within the ISM or amateur radio band. This applies to all ISM or amateur radio bands within the specified frequency range.			
e) The test may also be conducted at other baseband frequencies as determined by the RISK MANAGEMENT PROCESS.			
f) ME-DEVICES and ME-SYSTEMS with a DC power input, designed for operation with a transformer with an AC input and a DC output, must be tested with a transformer corresponding to the parameters specified by the MANUFACTURER of the ME-DEVICE or ME-SYSTEM. IMMUNITY TEST LEVELS must be applied to the transformer AC input.			
g) Applies only to ME-DEVICES and ME-SYSTEMS designed for connection to single-phase AC mains.			
h) Slash-separated data, for example 10/12, means 10 cycles at 50 Hz or 12 cycles at 60 Hz line frequency.			

- i) ME-DEVICES and ME-SYSTEMS with a RATED input current of more than 16 A per phase require a single interruption for 250/300 cycles at any angle and for all phases simultaneously (if applicable). ME-DEVICES and ME-SYSTEMS with a buffer mode after testing must return to the mains-powered operation. ME-DEVICES and ME-SYSTEMS with a RATED input current less than or equal to 16 A require interruption on all phases.
- j) ME-DEVICES and ME-SYSTEMS without transitional protective devices in the primary mains circuit should be tested only with ± 2 kV between the power supply wires and the earth and with ± 1 kV between individual power supply devices.
Does not apply to ME-DEVICES and ME-SYSTEMS with PROTECTION CLASS II.
It is necessary to use a direct connection.
- m) Root-mean-square value before modulation is applied.
- n) The ISM bands (frequency ranges used for industrial, scientific and medical purposes) between 0.15 MHz and 80 MHz are 6.765 MHz - 6.795 MHz; 13.553 MHz - 13.567 MHz; 26.957 MHz - 27.283 MHz; and 40.66 MHz - 40.70 MHz. The amateur radio bands between 0.15 MHz and 80 MHz are 1.8 MHz - 2.0 MHz, 3.5 MHz - 4.0 MHz, 5.3 MHz - 5.4 MHz, 7 MHz - 7.3 MHz, 10.1 MHz - 10.15 MHz, 14 MHz - 14.2 MHz, 18.07 MHz - 18.17 MHz, 21.0 MHz - 21.4 MHz, 24.89 MHz - 24.99 MHz, 28.0 MHz - 29.7 MHz and 50 MHz - 54.0 MHz.
- o) Applicable for ME-DEVICES and ME-SYSTEMS with RATED input current less than or equal to 16 A per phase, and for ME-DEVICES and ME-SYSTEMS with RATED input current more than 16 A per phase.
- p) Applicable for ME-DEVICES and ME-SYSTEMS with RATED input current less than or equal to 16 A per phase.
- q) For ME-DEVICES equipped with a mains input transformer, the use of such a test at certain phase angles can lead to the activation and deactivation of the overcurrent protection device. This can happen due to magnetic saturation of the transformer core after a voltage dip. If this happens, it is necessary to ensure the BASIC SAFETY of the ME-DEVICE or ME-SYSTEM during and after the test.
- r) For ME-DEVICES and ME-SYSTEMS, which have a large number of possible mains voltage settings or are equipped with an automatic mains voltage selection system, tests should be carried out with the mains voltage values corresponding to the minimum and maximum RATED VALUE of the mains input voltage, which is less than 25% of the maximum mains input voltage within the RATED VALUE range, the test must be carried out with the mains input voltage within this range. For corresponding calculation examples see Table 1 Note c).

Table 4: The IMMUNITY and SHELL test provisions for HF communications devices (according to IEC60601-1-2).

Test frequency (MHz)	Frequency range a) (MHz)	Radio communication ^{a)}	Modulation ^{b)}	Maximum power (W)	Distance (m)	TEST LEVEL OF IMMUNITY (V/m)
385	380 - 390	TETRA 400	Pulse modulation b) 18 Hz	1.8	0.3	27
450	430 - 470	GMRS 460, FRS 460	FM ^{c)} ± 5 kHz step 1 kHz sine	2	0.3	28
710	704 - 787	LTE band 13, 17	Pulse modulation b) 217 Hz	0.2	0.3	9
745						
780						
810	800 - 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE band 5	Pulse modulation b) 18 Hz	2	0.3	28
870						
930						
1720	1700 - 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE band 1, 3, 4, 25; UMTS	Pulse modulation b) 217 Hz	2	0.3	28
1845						
1970						
2450	2400 - 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE band 7	Pulse modulation b) 217 Hz	2	0.3	28
5240	5100 - 5800	WLAN 802.11 a/n	Pulse modulation b) 217 Hz	0.2	0.3	9
5500						
5785						
NOTE: If necessary, to achieve the TEST LEVEL OF IMMUNITY, the distance between the transmitting antenna and the ME-DEVICE or ME-SYSTEM can be reduced to 1 m. A test distance of 1 m is permitted in accordance with IEC 61000-4-3.						
a) For some radio services, only the radio frequencies for the mobile communications device with the 8en: uplink base station are included in the table.						
b) The carrier frequency should be modulated with a square wave signal with a duty cycle of 50%.						
c) As an alternative to frequency modulation (FM), pulse modulation with a duty cycle of 50% from 18 Hz can be used, since this, although it is also not an actual modulation, still represents the worst case.						

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KIRCHNER & WILHELM GmbH + Co. KG