



Martillos de reflejo

Reflex Hammers



COLORFLEX



147

Colorflex

– martillo de reflejo

– grande

– negro

– con punta palpadora neurológica

– reflex hammer

– large

– black

– with neurological sensor tip

aprox. | approx. 22 cm

aprox. | approx. 100 g

REF 07.15204.021

Colorflex

– martillo de reflejo

– grande

– azul

– con punta palpadora neurológica

– reflex hammer

– large

– blue

– with neurological sensor tip

aprox. | approx. 22 cm

aprox. | approx. 100 g

REF 07.15204.031

Colorflex

– martillo de reflejo

– grande

– gris

– con punta palpadora neurológica

– reflex hammer

– large

– grey

– with neurological sensor tip


aprox. | approx. 22 cm

aprox. | approx. 100 g

REF 07.15204.262

BABINSKY

ROSSIER

**Babinsky**

- martillo de reflejo
- con mango metálico, intercambiable por empalme roscado

- reflex hammer
- metal handle, interchangeable with screw fitting

Ø 45 mm
aprox. | approx. 21 cm
aprox. | approx. 85 g

REF 07.05100.001

Babinsky

- martillo de reflejo
- con mango elástico de material sintético, intercambiable por empalme roscado

- reflex hammer
- elastic plastic handle interchangeable with screw fitting

Ø 45 mm
aprox. | approx. 23 cm
aprox. | approx. 70 g

REF 07.05300.001

Rossier

- martillo de reflejo
- con mango elástico de material sintético

- reflex hammer
- elastic plastic handle

Ø 50 mm
aprox. | approx. 24 cm
aprox. | approx. 68 g

REF 07.10324.001

Rossier

- martillo de reflejo
- con mango elástico de material sintético

- reflex hammer
- elastic plastic handle

Ø 50 mm
aprox. | approx. 35 cm
aprox. | approx. 71 g

REF 07.10336.001

VARIOFLEX

**Varioflex**

- martillo combinado según Fassbender
- con mango metálico macizo con graduación en cm
- con pincel
- con aguja picanervios

- Fassbender combination hammer
- with solid metal handle, with scale setting in cm
- with brush
- with needle

aprox. | approx. 22 cm

aprox. | approx. 150 g

REF 07.25130.001



**Trömner**

- martillo de reflejo
- con mango metálico, modelo original
- pesado
- sin pincel
- sin aguja

- reflex hammer
- with metal handle: original model
- heavy
- without brush
- without needle

aprox. | approx. 23 cm
aprox. | approx. 210 g

Trömner

- martillo de reflejo
- con mango de material sintético negro
- pesado
- con pincel
- sin aguja

- reflex hammer
- with black plastic handle
- heavy
- with brush
- without needle

aprox. | approx. 23 cm
aprox. | approx. 170 g

Trömner

- martillo de reflejo
- con mango de material sintético negro
- ligero
- con pincel
- con aguja

- reflex hammer
- with black plastic handle
- light
- with brush
- with needle

aprox. | approx. 23 cm
aprox. | approx. 160 g

Trömner

- martillo de reflejo
- con mango de material sintético negro
- pesado
- con pincel
- con aguja

- reflex hammer
- with black plastic handle
- heavy
- with brush
- with needle

aprox. | approx. 23 cm
aprox. | approx. 200 g

Trömner

- martillo de reflejo
- con mango de material sintético negro
- ligero
- con pincel
- con aguja

- reflex hammer
- with black plastic handle
- light
- with brush
- with needle

aprox. | approx. 23 cm
aprox. | approx. 130 g

WITROE



BUCK



Witroe

- martillo de reflejo
- con mango de material sintético
- pesado
- con pincel
- con aguja
- reflex hammer
- with black plastic handle
- heavy
- with brush
- with needle

aprox. | approx. 24 cm
aprox. | approx. 195 g

REF 07.40232.001

Buck

- martillo de reflejo
- con mango metálico
- ligero
- con pincel
- con aguja

- reflex hammer
- with metal handle
- light
- with brush
- with needle

aprox. | approx. 19 cm
aprox. | approx. 100 g

REF 07.45130.001

Buck

- martillo de reflejo
- con mango de material sintético
- ligero
- con pincel
- con aguja

- reflex hammer
- with plastic handle
- light
- with brush
- with needle

aprox. | approx. 19 cm
aprox. | approx. 80 g

REF 07.45230.001